

LABORHITS

Industrial Jerrycans, HDPE

With UN approval*



Capacity L	Colour	Width mm	Depth mm	Height mm	Weight g	Thread external diam. mm	PK	Cat. No.
2.0	natural	145	109	196	110	45	1	9.140 021
2.0	natural	145	109	196	110	45	25	4.679 547
2.5	natural	153	115	208	125	45*	1	9.140 022
2.5	natural	153	115	208	125	45*	25	4.679 548
3.0	natural	153	115	240	140	45*	1	9.140 023
3.0	natural	153	115	240	140	45*	25	4.679 550
5.0	natural	190	145	251	250	50*	1	9.140 025
10.0	natural	230	190	310	450	50*	10	4.679 546
5.0	natural	190	145	251	250	50*	16	4.679 551
6.0	natural	192	145	294	265	45	1	9.140 026
8.0	natural	230	190	270	320	45	1	9.140 027
10.0	natural	230	190	310	450	50*	1	9.140 030
20.0	natural	290	246	385	950	60*	1	9.140 037
20.0	natural	290	246	385	950	60*	6	4.679 549
20.0	blue	290	246	385	950	60*	1	9.140 032
25.0	blue	290	246	450	1100	60*	1	9.140 031
30.0	blue	320	284	435	1250	60*	1	9.140 033
30.0	natural	320	284	435	1250	60*	1	9.140 038

Matching screw caps, pump dispenser and wrenches on the back side.

LABORHITS

Accessories for Industrial Jerrycans, HDPE

Screw caps

Pump dispenser

Wrenches



For	Description	PK	Cat. No.
5 / 10 litre jerrycans	Black, with tear-off ring and vent	1	9.140 068
2 / 2.5 / 3 / 6 / 8 litre jerrycans	Black, with tear-off ring and vent	1	9.140 048
20 / 25 / 30 litre jerrycans	Black, with tear-off ring and vent	1	9.140 050
2 / 2.5 / 3 / 6 / 8 litre jerrycans	Black, with tear-off ring	1	9.140 051
20 / 25 / 30 litre jerrycans	Black, with tear-off ring	1	9.140 047
5 / 10 litre jerrycans	Black, with tear-off ring	1	9.140 046
2 / 2.5 / 3 / 6 / 8 litre jerrycans	Red, with tear-off ring	1	9.140 040
5 / 10 litre jerrycans	Red, with tear-off ring	1	9.140 042
20 / 25 / 30 litre jerrycans	Red, with tear-off ring	1	9.140 043
2 / 2.5 / 3 / 6 / 8 litre jerrycans	With stopcock	1	9.140 062
5 / 10 litre jerrycans	With stopcock	1	9.140 061
20 / 25 / 30 litre jerrycans	With stopcock	1	9.140 063
2 / 2.5 / 3 / 6 / 8 litre jerrycans	With pump dispenser	1	9.140 054
10 litre jerrycans	With pump dispenser	1	9.140 065
20 litre jerrycans	With pump dispenser	1	9.140 070
25 litre jerrycans	With pump dispenser	1	9.140 071
30 litre jerrycans	With pump dispenser	1	9.140 072
2 / 2.5 / 3 / 6 / 8 litre jerrycans	Wrench, HDPE to apply / remove	1	9.140 073
5 / 10 / 20 / 25 / 30 litre jerrycans	Wrench, HDPE to apply / remove	1	9.140 108
2 / 2.5 / 3 / 6 / 8 litre jerrycans	Red, with tear-off ring	10	9.140 012
5 / 10 litre jerrycans	Red, with tear-off ring	10	9.140 013
20 / 25 / 30 litre jerrycans	Red, with tear-off ring	10	9.140 014
5 / 10 litre jerrycans	Black, with tear-off ring	10	9.140 015